


Loading Guide

for MediceI ACCUJECT 2.2* with enVista®TORIC (MX60T)



Viscoelastic colored for visualization.

Entering from the loading chamber side, apply Amvisc® Plus, or other Bausch + Lomb viscoelastic, directly into the conical tip and then apply 2 thin lines into the lateral grooves within the loading chamber.



Advance the plunger tip to the outer edge of the cartridge wings as shown.




Open the vial containing the IOL and, using nonserrated forceps, remove the lens by grasping and carefully pulling it out vertically from the center slot at the top portion of the vial.



Lens colored for visualization.

Position the lens in the middle of the loading chamber so that the anterior side is up and the lens is in a reverse-S orientation.




Lens colored for visualization.

Apply slight downward pressure with the forceps to push the lens and haptics down to ensure they are properly seated under the grooves.




Lens colored for visualization.

Slightly close the cartridge wings to hold the lens in place and then advance the plunger so that the haptics are compressed. The compression is correct when the haptic is pointing toward, but not touching, the optic.



Next, close the cartridge wings together until the click-lock mechanism engages.



Push the lens into the conical tip by advancing the plunger until the spring just contacts the outer edge of the cartridge wings.



Pull the plunger back sufficiently to visually confirm that the lens remains in the conical tip and then push forward again. The lens is now loaded and ready for injection.



With the conical tip bevel facing down, inject the lens by applying continuous pressure on the plunger until the lens is fully expressed from the tip. Clockwise injector rotation will compensate for any lens rotation.

*MediceI ACCUJECT 2.2 1P Rev. 2

enVista®TORIC (MX60T) is not approved for use in the United States.

Precautions and Warnings: Before loading the injector, consult the enVista®TORIC Directions for Use (DFU) as well as the ACCUJECT 2.2 1P (Rev. 2) Directions for Use for further safety information. Do not resterilize or reuse the ACCUJECT 2.2 1P (Rev. 2) injector. Single use only. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Store the ACCUJECT injector between 5° and 40°C (41° and 104°F). Folding and compression of the lens should be done just prior to insertion and delivery. Failure to follow the instructions for use may result in patient injury.

©2012 Bausch & Lomb Incorporated. ®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates. Other brands/products are trademarks of their respective owners. SU6751 08/12

new enVista®TORIC
Hydrophobic Acrylic Intraocular Lens
the new measure in toric stability

BAUSCH + LOMB
See better. Live better.