



## Wired Power Supplies

Whether your power supply needs to be stationary or portable Keeler Wired Power Supplies are robust and reliable.

### SmartPack

SmartPack is a stylish power supply system that combines 3 options in one product. No longer will you need to purchase a portable system for domiciliary use, another system for mounting on your wall and another for tabletop use.

### Fast charge and long life

Nickel Metal Hydride technology provides a genuine long life and a fast charge. Instruments can be used when battery is charging. LEDs indicate when battery life is low, charging or fully charged.



SmartPack

### WallPack

Alternatively, if your power supply needs to be stationary you can purchase the WallPack which fits neatly on the wall. Your indirect can be mounted on the bracket it comes with.

Which system is right for you?

Vantage Plus Usable Life		
Power Source	Natural LED	Bulb
Slimline Battery	6 hours*	1.5 hours*
Standard Battery	11 hours*	2.5 hours*
SmartPack	17 hours*	4.25 hours*

\* This measurement is based upon the illumination being used at normal examination power

### Part numbers

1945-P-1001	SmartPack
EP39-22079	SmartPack Battery
1945-P-1000	WallPack



WallPack

For Volk Indirect Lenses  
see pages 61-65



*Slimline Lithium Polymer Battery*



*Slimline Wireless Charger*

## Wireless Power Supplies

Keeler wireless technology allows you to move from room to room without any constraints. The chargers are either desktop or wall mountable and store an extra lithium battery as a back up, should you need it.

### New Slimline Lithium Polymer Battery

The New Slimline Lithium Polymer battery can give you up to 6 hours of use, weighs only 53 grams, and is nearly 70% lighter than the standard battery. It is also much smaller and can fit nicely into the palm of your hand.

### New Slimline Wireless Charger

The new wall mountable charger remains as durable as the Standard Wireless battery charger with the facility to store your extra battery in the top compartment. The charger has an easy release button to remove your extra battery.

### Standard Lithium Battery

Our original lithium battery lasts up to 11 hours of use and fits neatly into the charging unit.



### Standard Wireless Charger

The Standard Wireless Charger can be wall or desk mounted with space for a spare battery, which means you will always have a back up. An LED indicates when the battery needs charging.



*Slimline Lithium Polymer Battery*



*Standard Lithium Battery*

### Part numbers

1919-P-5338 Slimline Lithium Polymer Battery

1945-P-5012 Slimline Wireless Charger

1919-P-1013 Standard Lithium Battery

1941-P-1384 Standard Wireless Charger